

Endocavity

Ultrasound Scanning & Biopsy

Probe Covers



Developers & Manufacturers

Private Labeling

Custom Packaging

Certified Cleanrooms

Sterile & Non-sterile

Protek's probe covers are a proven shield against cross contamination and protects your equipment from blood and other contaminants from getting into hard to clean areas. Tested to isolate a site of surgical incision from microbial or other contaminants.

Soft Non-latex:

- Recommended alternative to latex
- Soft and comfortable material, emulates the properties of human tissue for distortion-free scanning
- Stretchable and available in many sizes

ORDER INFORMATION

Soft Non-latex

2201	4"-1.25"x12"	(10.2-3.2x30.5cm)	Non-sterile cover and bands (no gel), 50/Box
2220	5"-2"x24"	(12.7-5.1x61cm)	Sterile cover, bands and sterile gel, 24/Box
2225	5"-2"x24"	(12.7-5.1x61cm)	Non-sterile cover and bands (no gel), 24/Box

Latex

3230	1.4"x7.9"	(3.5x20cm)	Sterile cover with bands, 40/Box
3231	1"x7.9"	(2.6x20cm)	Non-sterile covers with bands, 24/Box
3232	0.8"x7.9"	(2.0x20cm)	Non-sterile covers with bands, 24/Box
3233	0.8"x7.9"	(2.0x20cm)	Sterile covers with bands, 24/Box
3235	1.4"x7.9"	(3.5x 20cm)	Non-sterile cover with bands, 40/Box
3238	1.4"x7.9"	(3.5x 20cm)	Non-sterile covers bulk, 100/Box



For nearest distributor & order information contact:
4125 Westcor Court, Coralville, IA 52241, USA
Phone: (319)545-7100 • Fax: (319)545-7099
E-mail: info@protekmedical.com
Web: www.protekmedical.com