

General Purpose

Ultrasound Scanning & Biopsy

Probe Covers



Developers & Manufacturers

Private Labeling

Custom Packaging

Certified Cleanrooms

Sterile & Non-sterile

Protek's probe covers are a proven shield against cross contamination and protects your equipment from blood and other contaminants from getting into hard to clean areas. Tested to isolate a site of surgical incision from microbial or other contaminants.

Soft & Clear Non-latex:

- Recommended alternative to latex
- Soft and comfortable material, emulates the properties of human tissue for distortion-free scanning
- Available in many sizes

ORDER INFORMATION

Soft Non-latex Sterile Probe Cover Kits (24/box)

2097	5.5"x48"	(14x121.9cm)	3D end with <i>PullUp</i> ™, sterile gel and bands
2099	5.5"x24"	(14x61cm)	3D end with <i>PullUp</i> ™, sterile gel and bands
2230	3.5"x24"	(8.9x61cm)	With sterile gel and bands
2250	5.5"x24"	(14x61cm)	3D end with sterile gel and bands
2270	7"x24"	(17.8x61cm)	With sterile gel and bands
2320	5.5"x36"	(14x91.4cm)	3D end with sterile gel and bands
2450	5.5"x48"	(14x121.9cm)	3D end with sterile gel and bands

Clear Non-latex Sterile Probe Cover Kits (24/box)

1096	5"x48"	(12.7x121.9cm)	With <i>PullUp</i> ™, sterile gel and bands
1240	4"x24"	(10.2x61cm)	With sterile gel and bands
1460	6"x48"	(15.2x121.9cm)	With sterile gel and bands

PullUp™ is a trademark & patented by Protek Medical Products, Inc.



For nearest distributor & order information contact:
4125 Westcor Court, Coralville, IA 52241, USA
Phone: (319)545-7100 • Fax: (319)545-7099
E-mail: info@protekmedical.com
Web: www.protekmedical.com