


ProV-Access™

Ultrasound Biopsy
Needle Guides
For Vascular Access Procedures



U
L
T
R
A
S
O
U
N
D
A
C
C
E
S
S
O
R
I
E
S

Developers & Manufacturers

Private Labeling

Custom Packaging

Certified Cleanrooms

Sterile & Non-sterile

A vital physician's tool to be used as an accessory to ultrasound transducers, guiding a needle or catheter into a target to perform a biopsy or precise needle placement.

Sterile • Disposable • Non-latex

- For vascular access ultrasound guided procedures
- Fits multiple gauges of needles (18 through 21)
- Patented needle release for tilt out procedures
- Clear pad to visualize needle entry
- Easily secures to the probe

ORDER INFORMATION

Contains: needle guides, probe cover, sterile gel and bands
Packaged with the *PullUp™* Probe Cover:
easy to open • easy to add gel • easy to pull up over the probe

3 Depth Sterile Procedure Kits 10/Box

#5410	1.0cm, 1.5cm, 3.5cm
#5414	1.0cm, 1.5cm, 2.0cm
#5415	1.5cm, 2.0cm, 3.5cm

Single Depth Sterile Needle Guides 10/Box

#5404	0.5cm
#5413	1.0cm
#5412	1.5cm
#5405	2.0cm
#5411	3.5cm

One or more U.S. patents apply



For nearest distributor & order information contact:
4125 Westcor Court, Coralville, IA 52241, USA
Phone: (319)545-7100 • Fax: (319)545-7099
E-mail: info@protekmedical.com
Web: www.protekmedical.com

PROTEK