

# RetroFit™

Ultrasound Biopsy  
Needle Guides



U  
L  
T  
R  
A  
S  
O  
U  
N  
D  
A  
C  
C  
E  
S  
S  
O  
R  
I  
E  
S

Developers & Manufacturers

Private Labeling

Custom Packaging

Certified Cleanrooms

Sterile & Non-sterile

*A vital physician's tool to be used as an accessory to ultrasound transducers, guiding a needle or catheter into a target to perform a biopsy or precise needle placement.*

### Sterile • Disposable • Non-latex

- Attaches securely to mounting bracket
- Color coded to match needle gauge
- Clear for visual needle placement
- Easy quick release pad for needle
- Needle gauge specific
- No assembly required
- 13 thru 23 gauge

## ORDER INFORMATION

Contains: needle guide, soft probe cover, sterile gel and bands

### Retrofit™ Sterile Needle Guide Probe Cover Kits 24/Box

Gauge	13	14	15	16	17	18	19	20	21	22	23
Order#	5513	5514	5515	5516	5517	5518	5519	5520	5521	5522	5523

One or more U.S. patents apply  
Packaged with the PullUp™ Probe Cover:  
easy to open • easy to add gel • easy to pull up over the probe



For nearest distributor & order information contact:  
4125 Westcor Court, Coralville, IA 52241, USA  
Phone: (319)545-7100 • Fax: (319)545-7099  
E-mail: [info@protekmedical.com](mailto:info@protekmedical.com)  
Web: [www.protekmedical.com](http://www.protekmedical.com)